

Anti-Rat IgA (MARA-1)

Mouse Monoclonal Antibody

Cat. No.	Conjugate	Unit Size
MARA-1-P-0.1	purified	100µl
MARA-1-P-0.5	purified	500µl
MARA-1-P-1	purified	1ml
MARA-1-HRP-0.5	Peroxidase	500µl
MARA-1-HRP-1	Peroxidase	1ml
MARA-1-FITC-0.5	FITC	500µl
MARA-1-FITC-1	FITC	1ml
MARA-1-BIOT-0.5	Biotin	500µl
MARA-1-BIOT-1	Biotin	1ml

Product Information

Supplied as: 1 mg/ml affinity purified antibody prepared in PBS buffer, pH 7.6, with 0.09% sodium azide as preservative FITC and Biotin labelled antibodies are supplied in PBS buffer, pH 7.6, with 0.09% sodium azide and 50% glycerol.

Origin: Hybridoma derived by hybridization of spleen cells from immunized mice with the Sp2/O mouse myeloma cell line.

Immunogen: IgA (IR1060, IR22, IR699) (1)

Clone: MARA-1

Host / Isotype: Mousekappa IgG1

Reactivity: Rat alpha heavy chain of immunoglobulin (determined by immunodot). Recognizes a conformational epitope in Western Blott. Does not recognize the secretory piece.

Cross-Reactivity: not determined

Avidity: on rat IgA: $1.1 \times 10^{10} \text{M}^{-1}$

Applications:

- CAN BE USED AS SECOND ANTIBODY IN IMMUNOASSAYS
- CAN BE USED FOR CAPTURE ELISA: GOOD BINDING ON PLASTICS
- CAN BE USED IN WESTERN BLOT ANALYSIS (only in non denaturated conditions)
- CAN BE USED IN IMMUNOAFFINITY CHROMATOGRAPHY (Solid phase Sepharose 4B CNBr act.)
- CAN BE COATED ON NITROCELLULOSE (DOT-ELISPOT)

Storage and Stability: For longterm storage store product at -20°C.

At 2-8°C the undiluted stock solution of the antibody is stable for a few weeks. Repeated freezing and thawing of unconjugated antibody should be avoided. It is recommended to aliquot unconjugated antibody before freezing. Labelled antibodies (Peroxidase, FITC, Biotin) contain 50% glycerol as antifreezing and remain liquid at -20°C. In this case aliquoting is not required.

It's recommended to prepare a working dilution fresh prior to use.

Expiry information is given on the label and / or data sheet supplied upon delivery

Anti-Rat IgA (MARA-1)

Mouse Monoclonal Antibody

Safety Information: Standard Laboratory Practices should be followed when handling this material. The chemical, physical, and toxicological properties of this material have not been thoroughly investigated. Purified, FITC, and Biotin labelled antibody contains 0.09% sodium azide (a poisonous and hazardous substance) as preservative. Although this concentration is not regarded as dangerous to health, appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. Peroxidase conjugate contains no hazardous constituents.

Note: This information sheet contains general product information only. For detailed lot specific information consult the label and data sheet supplied by the manufacturer upon delivery of the product.

References: (1) Bazin H., Beckers A. & Querinjean P. Three classes and four subclasses of rat immunoglobulins: IgM, IgA, IgE and IgG1, IgG2a, IgG2b, IgG2c. Eur. J. Immunol. 1974,4:44-48.
Selected references of scientific publications in which MARA-1 MAb was used
Soares et al. Transplantation, 1993, 56: 1427-1433
Vanhove B., Bazin H. Mol. Immunology, 1992, 29:1-8
Vanhove B. & Bazin H. Immunology 1993, 79: 580-586.
Latinne D., et al. Immunol. Rev. 1994, 141: 95-125
Soares et al., Transplant. Proc. 1995, 27: 282-285

For more information, see:

"Rat Hybridomas and Rat Monoclonal Antibodies". H. BAZIN (Ed.). CRC Press, Boca Raton, Florida, USA, 1990, 515 pages.

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.

This product is warranted to perform in conformance with the product specifications provided by the manufacturer and to be free from defects in material and workmanship. Products are supplied for research use only. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated, this warranty is limited to the period of stability as indicated on the label and / or data sheet at the time of delivery of the product when the product is subjected to normal, proper and intended usage. The remedy for non-conforming products is limited to replacement of or refund for the non-conforming product(s). There is no obligation to replace products as the result of accident or inappropriate use, improper storage or handling of the products.

Manufacturer:

SynAbs SA
6 Chemin du Cyclotron
B-1348 Louvain-la-Neuve, Belgium

Distribution and Support:

dianova GmbH
Warburgstr. 45
20354 Hamburg
www.dianova.com

Technical Support:
Phone: 040 45 06 70
Email: info@dianova.de

Ordering Information:
Phone: 040 45 06 70
Fax: 040 45 06 74 90
Email: order@dianova.de